

PMT \_\_\_\_\_ HV board \_\_\_\_\_ HV module \_\_\_\_\_

Data directory name \_\_\_\_\_

Your name \_\_\_\_\_

Gain calibration Date \_\_\_\_\_ Time \_\_\_\_\_ Room temperature \_\_\_\_\_

### Settings

- Helmholtz coil X 18.0 V Y 12.0 V
- Low voltage \_\_\_\_\_
- Current when  $V_{ctrl} = 3.5$  V \_\_\_\_\_
- Filter 1 % + 0.1 %
- Intensity 5.8
- Frequency 100 Hz
- Resolution 500 ps
- Horizontal range (MUST include  $-100$  ns  $- \sim 150$  ns) \_\_\_\_\_
- Trigger low enough
- Scan control voltages ( ) ( ) ( ) ( )

### Analyses

Gain scan

- Gate start time \_\_\_\_\_
- $V_{ctrl}$  @1e7 gain \_\_\_\_\_

$V_{ctrl}$  @1e7 gain

- SPE mean charge \_\_\_\_\_
- SPE peak height \_\_\_\_\_

Linearity & afterpulse Date \_\_\_\_\_ Time \_\_\_\_\_ Room temperature \_\_\_\_\_

### Settings

- Helmholtz coil X 18.0 V Y 12.0 V
- Low voltage \_\_\_\_\_
- Control voltage (MUST  $V_{ctrl}$  @1e7 gain) \_\_\_\_\_
- Current \_\_\_\_\_
- Intensity 14
- Frequency 10 Hz
- Resolution 1 ns
- Horizontal range  $-20 \mu s - 20 \mu s$

### Analyses

2D scan Date \_\_\_\_\_ Time \_\_\_\_\_ Room temperature \_\_\_\_\_

### Settings

- Helmholtz coil X 18.0 V Y 12.0 V
- Low voltage \_\_\_\_\_
- Control voltage (MUST  $V_{ctrl}$  @1e7 gain) \_\_\_\_\_
- Current \_\_\_\_\_
- Filter 1 % + 1 %
- Intensity 8.5 (Pulse height should  $\sim$  600 mV)
- Frequency 100 Hz
- Resolution 500 ps
- Horizontal range (MUST include  $-100$  ns  $- \sim$  150 ns) \_\_\_\_\_
- #waveforms at each point \_\_\_\_\_

### Analyses

■

### Dark rate

Put into freezer Date \_\_\_\_\_ Time \_\_\_\_\_

Start taking data Date \_\_\_\_\_ Time \_\_\_\_\_ Freezer temperature (MUST  $\sim -30$  °C) \_\_\_\_\_

Take out from freezer Date \_\_\_\_\_ Time \_\_\_\_\_ Freezer temperature \_\_\_\_\_

### Settings

- Black tape Yes No
- Low voltage \_\_\_\_\_
- Control voltage (MUST  $V_{ctrl}$  @1e7 gain) \_\_\_\_\_
- Current \_\_\_\_\_
- Resolution 2 ns
- Horizontal range  $-20$  ms  $- 20$  ms
- Vertical range 2 mV/div
- Trigger 2 mV

### Analyses

- #files (MUST  $\geq 500$ ) \_\_\_\_\_
- Gaussian fitting range \_\_\_\_\_
- SPE mean charge \_\_\_\_\_
- Dark rate for 0.25 PE threshold and 100 ns deadtime \_\_\_\_\_

Put the json file into database

Print out the plots